

Excellent Surgical Instruments

EXCELLENT QUALITY FOR O.P
DISPOSABLE QUALITY FOR FLOOR GRADE



Excellent Instruments Are Manufactured In Pakistan With Modern Technical Expertise, From The Finest Quality Surgical Steel Which Is Heat Treated To Suit The Quality Needed For The Various Surgical Instruments

Director By : Malik Shahid Saleem

Our Instruments Are Finished To Ensure The Highest Degree Of Resistance To Rust And Stains. The Modern Technique Used Produces An Instruments Of Both Quality And Durability, Which Will Stand Up To Normal Use. They Have Been Found To Be Economical For C.S.R And General Medical Practice.



www.salim-group.com

Dear Customer :

We are happy to present the ninth edition of our catalogue. We have improve on it in order to make it Simpler and easier to follow. The selection of instruments has also been increased. This we have done because of the quality workmanship put into the manufacturer of Excellent surgical instruments, and because of the wide acceptance of these instruments in the hospitals. In many cases our surgical instruments can replace other well known brands at remarkable savings.

With the modern trend of computerization, we feel that the time has come to etch each instrument with its catalogue number. This will help inventory control tremendously, and it will also make it easier for hospital staff to order and sort the instruments.

All rights reserved. This catalogue or any part thereof, shall not be reproduced in any form without written permission from Salim Group Of Industries.

First Edition.


January 2013

INDEX


Page/Column:	Item:	Item:	Page/Column:
Alligator Forceps	6/2	Nail Nippers	7/3
Angiotribes	6/2	Nail Scissors	8/1-8/2
Applying And Removing Forceps	4/1	Nasal Specula	6/2
Artery Clamps	5/1	Needle Holders	5/1-3
Bandage Scissors	1/3	Operating Scissors	2/1-3
Babinsky Neurological Hammer	1/1	Percussion Hammers	1/1
Bone Cutting Forceps	7/1	Plaster Shears/knives	2/1
Bone Sawa	8/3	Probes	7/1
Buck Neurological Hammers	1/1	Retractors	6/3
Chiropody Instruments	7/3-8/1	Scalpel Handles	7/2
Comedone Extractors	1/1	Scissors	2/1-4
Corn Cutter With Blade	8/1	Splinter Forceps	3/3
Dental Instruments	7/3	Sponge Forceps	6/1
Diagnostic Instruments	1/1	Sterilizing Instruments	1/2-3
Dieffenbah Serrefines	5/1	Stitch Scissors	2/3
Dix Needle And Spud	1/1	Suction Tubes	8/3
Dressing Forceps	3/2	Surgical Blades	9/3
Ear Instruments	6/2	Suture And Wire Scissors	2/3
Eye Dressing Forceps	3/2	Swiss Style Forceps	3/3-4/1
Finger Nail Drill	1/1	Tenaculum Forceps	6/1
Finger Ring Cutter	1/2	Tissue Forceps	3/2
Forceps	3/2	Toenial Splitter	7/3
General Purpose Scissors	9/1	Tongue Forceps	1/2
Gynecological Instruments	6/1	Tooth Extractor	7/3
Grooved Direcors	7/2	Towel Clamps	5/3
Hemostatic Forceps	4/2	Tube Gauze Applicators	9/3
Intestinal Forceps	4/1	Tube Occluding Forceps	8/3
Introducing Forceps	1/2	Uterine Dilators	6/1
Iris Forceps	3/3	Uterine Dressing Forceps	6/1
Micro Dissecting Scissors	3/1	Uterine Sounds	6/1
Miscellaneous Instruments	9/1	Utilty Scissors	1/3
Mouth Gags	6/3	Wartenberg Pinwheel	1/1

GUARANTEE - Our Surgical instruments are guaranteed to render satisfactory service provided they are used for the right purpose of operation with respect to their quality and design. Any instruments found defective in workmanship or material will be repaired or replaced at our discretion free of charge.


DIAGNOSTIC INSTRUMENTS




SGI-101 Taylor Percussion Hammer
18.5 cm, 7 1/4", Stainless
SGI-102 Taylor Percussion Hammer
18.5cm, 7 1/4", c.p.



SGI-103 Queens Percussion Hammer



SGI-104 Telescopic Percussion Hammer



SGI-105 Buck Neurological hammer
18cm, 7", stainless



SGI-106 Babinsky Neurological hammer,
25 cm, 10"



SGI-107 Wartenberg Pinwheel
stainless



SGI-108 Unna Comedone Extractor D/E



SGI-109 Comedone Extractor
With spear point




SGI-110 Schamberg Comedone Extractor
open end, stainless.




SGI-111 Schamberg Comedone Extractor
small end crimped, stainless




SGI-112 Ralks Finger nail Drill
with 3 drill points c.p




SGI-113 Dix Needle and Spud
reversible handle, stainless




SGI-114 Collins Tongue Seizing Forceps
small ring Jaws, 11.5cm, 4 1/2", s.s
SGI-115 Collins Tongue seizing forceps
large ring jaes, 18cm 7", s.s




SGI-116 Young Tongue
Seizing forceps
Soft Rubber jaws,
16.5 cm, 6 1/2", s.s




SGI-117 Magill Catheter introducing Forceps
24.5 cm, 9 3/4, Adult, Stainless
SGI-118 Magill Catheter introducing Forceps
20 cm, 8" Child, stainless



SGI-119 Finger Ring Cutter, stainless




SGI-120 Finger Ring Cutter, stainless
SGI-121 R/wplacement Blade for above




SGI-122 Finger Ring Cutter, (first aid Patted).

STERILIZING INSTRUMENTS



SGI-123 Sterilizing Forceps
utility straight 20 cm, 8", stainless
SGI-124 Sterilizing Forceps
utility curved 20 cm, 8", stainless
SGI-125 Sterilizing Forceps
utility straight 28 cm, 11", stainless
SGI-125 Sterilizing Forceps
utility curved 28 cm, 11", stainless



SGI-127 Sterilizing Forceps
three prong 20.5 cm, 8", stainless
Sterilizing forceps

SGI-128
three prong 30.5cm, 12, stainless



SGI-129 Cheatle Sterilizing Forceps
28 cm, 11", Stainless



SGI-30 Ralks Sterilizing forceps
for syringes 15 cm, 6" stainless
SGI-31 Ralks Sterilizing forceps
for syringes 25.5 cm, 10"

BANDAGE AND PLASTER INSTRUMENTS




SGI-132 Lister Bandage Scissors
9cm, 3 1/2", stainless

SGI-133 Lister Bandage Scissors
11.5cm, 4 1/2, stainless

SGI-134 Lister Bandage Scissors
14 cm, 5 1/2", stainless

SGI-135 Lister Bandage Scissors
18.5cm, 7 1/4, stainless



SGI-136 lister Bandage Scissor
one large ring 21.5cm, 8 1/4, Stainless




**SGI-137 KNOWLE BANDAGE SCISSORS
14CM, 5 1/2, STAINLESS**




SGI-138 Utility all Purpose Scissors
stainless{Available in all sizes}



SGI-139 first aid shears 21 . 6 cm , 8 1/2"s.s.



SGI-140 Utility All Purpos Nurses Scissors
14.5cm , s.s.



SGI-141 Bruns Plaster Shears
one blade serr. 27 cm ,10 3/4 ,
stainless



SGI-142 Stille Plaster Shears
26 cm ,10 1/2 " , stainless
SGI -142A Stille Plaster Shear
37 cm , 14 1/2, stainless



SGI-143 Esmarch Plaster Knife
18 cm , 7" , stainless




SGI-144 Bergmann Plaster Knife
18 cm , 7" , stainless



SGI-145 Reiner Plaster Knife
18 cm , 7" , stainless

Operating Scissor



SGI-146 Operating Scissors
straight, 11.5cm, 4 1/2 " B/B, s.s.
SGI-147 Operating Scissors
straight, 112.5cm, 5" , B/B, s.s.
SGI-148 Operating Scissors
straight, 14cm, 5 1/2 " , B/B, s.s.
SGI-149 Operating Scissors
straight, 15cm, 6 " , B/B, s.s.




SGI-150 Operating Scissors
straight, 18.5cm, 4 1/2 " S/S, s.s.

SGI-151 Operating Scissors
straight, 12.5cm, 5 " ,S/S, s.s.

SGI-152 Operating Scissors
straight, 14cm, 5 1/2 " S/S, s.s.

SGI-153 Operating Scissors
straight, 15cm,6 " S/S, s.s.




SGI-154 Operating Scissors
straight, 11.5cm,4 1/2 " S/S, s.s.


SGI-155 Operating Scissors
straight, 12.5cm,5 1/2 " S/S, s.s.

SGI-156 Operating Scissors
straight, 14cm,5 1/2 " S/S, s.s.


SGI-157 Operating Scissors
straight, 15cm,6 " S/S, s.s.




SGI-158 Enterotomy Scissors
21.5 cm, 8 1/2" , stianless




SGI-159 Iris Scissors
straight, 11.5cm,4 1/2 , s/s.




SGI-160 Iris Scissors
curved, 11.5cm,4 1/2, s/s.



SGI-161 Strabismus Scissors
straight, 12cm,4 3/4, s/s.



SGI-162 Strabismus Scissors
curved, 12cm,4 3/4, s/s.




SGI-163 Operating Scissors
curved, 11.5cm,4 1/2,B/B s/s.

SGI-164 Operating Scissors
curved, 12.5cm,5,B/B s/s.

SGI-165 Operating Scissors
curved, 14cm,5 1/2,B/B s/s.

SGI-166 Operating Scissors
curved, 15cm,6,B/B s/s.




SGI-167 Operating Scissors
curved, 11.5cm,4 1/2,S/B s/s.

SGI-168 Operating Scissors
curved, 12.5cm,5 " ,S/B s/s.

SGI-169 Operating Scissors
curved, 14cm,5 1/2,S/B s/s.

SGI-170 Operating Scissors
curved, 15cm,6,S/B s/s.




SGI-171 Operating Scissors
curved, 11.5cm,4 1/2,S/B s/s.

SGI-172 Operating Scissors
curved, 12.5cm,5 " ,S/B s/s.


SGI-173 Operating Scissors
curved, 14cm,5 1/2,S/B s/s.

SGI-174 Operating Scissors
curved, 15cm,6,S/B s/s.




SGI-175 Operating Stitch Scissors
Delicate,9 cm,3 1/2",Stainless

SGI-176 Operating Stitch Scissors
11.5 cm, 4 1/2" , stainless



SGI-177 Littauer junior Stitch Scissors
11.5 cm,4 1/2",Stainless

SGI-178 Littauer junior Stitch Scissors
14 cm,5 1/2",Stainless



SGI-179 Suture and wire Scissors
angular,12.5 cm,5",Stainless




SGI-180 Mayo Scissors
straight 14 cm,5 1/2",Stainless

SGI-181 Mayo Scissors
straight 17 cm,6 1/2",Stainless




SGI-182 Mayo Scissors
curved, 17 cm,6 3/4",Stainless


Operating Scissor




SGI-184 Metzenbaum Scissors
straight 15 cm ,6", stainless
SGI-185 Metzenbaum Scissors
straight 18 cm ,7", stainless




SGI-186 Metzenbaum Scissors
curved 15 cm ,6", stainless
SGI-187 Metzenbaum Scissors
curved 18 cm ,7", stainless




SGI-188 Dean Tonsil Scissors
angular, one blade serrated
17 cm, 6 3/4", stainless




SGI-189 Croen & Collar Scissors
straight 11.5 cm , 4 1/2,s.s.
SGI-190 Croen & Collar Scissors
curved 11.5 cm , 4 1/2,s.s.




SGI-191 Dissecting Scissors
angular 11.5 cm , 4 1/2,s.s.




Blunt points
SGI-192 Micro Dissecting Scissors
fine , spring handle
9 cm , 3 1/2, stainless



Very fine points
SGI-193 Micro Dissecting Scissors
12 cm , 4 3/4", stainless



SGI-193 Noyes
11 cm , 4 1/2", stainless



SGI-195 Dressing Forceps,11.5 cm , 4 1/2", s.s
SGI-196 Dressing Forceps, 12.5 cm , 5", s.s
SGI-197 Dressing Forceps,14 cm ,5 1/2", s.s
SGI-198 Dressing Forceps,15 cm , 6", s.s
SGI-199 Dressing Forceps,20 cm , 8", s.s
SGI-200 Dressing Forceps,25.5 cm , 10", s.s
SGI-201 Dressing Forceps,30.5 cm , 12", s.s



SGI-202 McIndoe Dressing Forceps
15.25 cm , 6", stainless



SGI-203 Semken Dressing Forceps
fine, 12.5 cm , 5", stainless



SGI-204 Adson Dressing Forceps
fine ,12 cm ,4 3/4", stainless



SGI-205 Eye Dressing Forceps
straight ,10 cm ,4", stainless




SGI-206 Eye Dressing Forceps
half curved ,10 cm ,4", stainless




SGI-207 Eye Dressing Forceps
full curved ,10 cm ,4", stainless

Tissue Forceps



SGI-208 Tissur Forceps
1x2 teeth,11.5 cm ,4 1/2", stainless
SGI-209 Tissur Forceps
1x2 teeth,12.5 cm ,5", stainless
SGI-210 Tissur Forceps
1x2 teeth,14 cm ,5 1/2", stainless
SGI-211 Tissur Forceps
1x2 teeth,15 cm ,6", stainless
SGI-212 Tissur Forceps
1x3 teeth,11.5 cm ,4 1/2", stainless
SGI-213 Tissur Forceps
1x3 teeth,12.5 cm ,5", stainless
SGI-214 Tissur Forceps
1x3 teeth,14 cm ,5 1/2", stainless



SGI-215 Semken Tissue Forceps
1x2 teeth,12 cm ,4 3/4", stainless



SGI-216 Gillies Tissue Forceps,1x2 teeth
15 cm 6".



SGI-217 Adson Tissue Forceps
1x2 teeth,fine 12 cm ,4 3/4", s.s




SGI-218 Brown Adson Tissue Forceps
side grasping, 12 cm ,4 3/4", s.s



SGI-219 Iris Forceps
1x2 teeth, straight, 10 cm ,4", s.s



SGI-220 Iris Forceps 1x2 teeth
half curved, 10 cm ,4", stainless




SGI-221 Iris Forceps 1x2 teeth
full curved, 10 cm ,4", stainless


Splinter forceps




SGI-222 Splinter Forceps
straight, 9 cm ,3 1/2", stainless




SGI-223 Splinter Forceps
straight, 10 cm ,4", stainless



SGI-224 Splinter Forceps
straight, 11.5 cm ,4 1/2", stainless




SGI-225 Feilchenfeld Splinter Forceps
7.5 cm ,3", stainless




SGI-226 Feilchenfeld Splinter Forceps
11.5 cm ,4 1/2", stainless




SGI-227 Carmlt Splinter Forceps
straight, 11.5 cm ,4 1/2", stainless




SGI-228 Carmlt Splinter Forceps
angular, 11.5 cm ,4 1/2", stainless



SGI-229 Ralks Splinter Forceps
15 cm, 6" stainless



SGI-230 #3 wiss Style Forceps
11 cm, 4 1/2" stainless



SGI-231 #4 wiss Style Forceps
11 cm, 4 1/2" stainless



SGI-132 #5 Swiss Style Forceps
11 cm , 4 1/2 , stainless




SGI-133 #6 Swiss Style Forceps
11 cm , 4 1/2 , stainless




SGI-134 #7 Swiss Style Forceps
11 cm , 4 1/2 , stainless

Applying And Removing Forceps




SGI-135 Michel cliip Applying and
Removing Forceps, Box joint
12.5 cm 5", stainless




SGI-136 Michel cliip Applying forceps
12.5 cm 5", stainless


Intestinal Forceps




SGI-237 Allis Tissue Forceps
4x5 teeth,15 cm ,6", stainless




SGI-238 Allis Tissue Forceps
5x6 teeth,18.5 cm ,6", stainless




SGI-239 Allis Tissue Forceps
5x6 teeth,18.5 cm ,7 1/4", stainless



SGI-240 Babcock intestinal Forceps
16 cm ,6 1/4",stainless



SGI-241 Babcock intestinal Forceps
20 cm ,8",stainless




SGI-242 Doyen intestinal Forceps
straight, 2 cm , 9", stainless




SGI-243 Doyen intestinal Forceps
straight, 23 cm , 9", stainless


HEMOSTATIC FORCEPS




SGI-244 Halstead Mosquito Forceps
straight,12.5 cm , 5", stainless




SGI-245 Halstead Mosquito Forceps
curved,12.5 cm , 5", stainless




SGI-246 Kelly Forceps
straight,14 cm , 5 1/2", stainless




SGI-247 Kelly Forceps
curved,14 cm , 5 1/2", stainless




SGI-248 Rankin Forceps
straight,16 cm , 6 1/4", stainless



SGI-249 Rankin Forceps
curved,16 cm , 6 1/4", stainless




SGI-250 Crile Forceps
straight,14 cm , 5 1/2", stainless




SGI-251 Crile Forceps
curved,14 cm , 5 1/2", stainless




SGI-252 Crile Forceps
straight,16 cm ,6 1/4", stainless




SGI-254 Hartman Mosquito Forceps
straight, 9 cm 3 1/2", stainless




SGI-255 Hartman Mosquito Forceps
curved, 9 cm 3 1/2", stainless




SGI-256 Providence Hospital Forceps
straight, 14 cm 5 1/2", stainless



SGI-257 Providence Hospital Forceps
curved, 14 cm 5 1/2", stainless




SGI-258 Rocheser- Pean Forceps
straight, 14 cm 5 1/2", stainless




SGI-259 Rocheser- Pean Forceps
straight, 16 cm 6 1/4", stainless




SGI-260 Rocheser- Pean Forceps
straight, 18.5 cm 7 1/4", stainless




SGI-261 Rocheser- Pean Forceps
curved, 14 cm 5 1/2", stainless




SGI-262 Rocheser- Pean Forceps
curved, 16 cm 6 1/4", stainless



SGI-263 Rocheser- Pean Forceps
curved, 18.5 cm 7 1/4", stainless



SGI-264 Rocheser- Ochsner Forceps
straight, 16 cm 6 1/4", stainless




SGI-265 Rocheser- Ochsner Forceps
straight, 18.5 cm 7 1/4", stainless




SGI-266 Rocheser- Ochsner Forceps
straight, 20 cm 8", stainless




SGI-267 Rocheser- Ochsner Forceps
curved, 16 cm 6 1/4", stainless




SGI-268 Rocheser- Ochsner Forceps
curved, 18.5 cm 7 1/4", stainless




SGI-269 Rocheser- Ochsner Forceps
curved, 20 cm 8", stainless



SGI-270 KocherForceps
straight, 14 cm , 5 1/2", stainless




SGI-271 KocherForceps
curved, 14 cm , 5 1/2", stainless




SGI-272 Rocheser- Carmalt Forceps straight, 16 cm 6 1/4", stainless
SGI-273 Rocheser- Carmalt Forceps curved, 16 cm 6 1/4", stainless

ANGIOTRIBES AND ARTERY CLAMPS



SGI-274 ferguson Angiotribe straight, 16.5 cm 6 1/2", stainless
SGI-275 ferguson Angiotribe curved, 16.5 cm 6 1/2", stainless
SGI-276 ferguson Angiotribe straight, 19 cm 7 1/2", stainless
SGI-277 ferguson Angiotribe curved, 19 cm 7 1/2", stainless




SGI-278 Dieffenbach Serrefine straight, 2.54 cm, 1", stainless
SGI-279 Dieffenbach Serrefine straight, 3.8 cm, 1 1/2", stainless
SGI-280 Dieffenbach Serrefine straight, 5.7 cm, 2 1/4", stainless




SGI-281 Dieffenbach Serrefine curved, 2.54 cm, 1", stainless
SGI-282 Dieffenbach Serrefine curved, 3.8 cm, 1 1/2", stainless
SGI-283 Dieffenbach Serrefine curved, 5.7 cm, 2 1/4", stainless


NEEDLE HOLDERS




SGI-284 Mayo- Hegar Needle Holder 12.5 cm 5", stainless
SGI-285 Mayo- Hegar Needle Holder 15 cm 6", stainless
SGI-286 Mayo- Hegar Needle Holder 18 cm 7", stainless
SGI-287 Mayo- Hegar Needle Holder 20 cm 8", stainless



SGI-288 Baumgartner Needle holder 14 cm 5 1/2", stainless




SGI-289 Crile-Wood Needle holder 14 cm 5 1/2", stainless




SGI-290 Olsen-Hegar Needle holder 14 cm 5 1/2", stainless
SGI-290 Olsen-Hegar Needle holder 16.5 cm 6 1/2", stainless




SGI-292 Gillies Needle holder combined with scissors 16 cm 6 1/4", stainless




SGI-293 Derf Needle holder 12 cm 4 3/4", stainless




SGI-294 Webster Needle holder smooth jaw 12.5 cm 5", s.s.



SGI-295 Halsey Needle holder smooth jaw 12.5 cm 5", s.s.




SGI-296 Mathieu Needle Holder for delicate sutures 14 cm , 5 1/2", s.s.



SGI-297 Mathieu Needle Holder with or without groove 14 cm , 5 1/2", s.s.

SGI-298 Mathieu Needle Holder with or without groove 17 cm , 6 3/4", s.s.




SGI-299 Mathieu Needle Holder with recessed ratchet 17 cm , 6 3/4", s.s.

SGI-300 Mathieu Needle Holder with recessed ratchet 20 cm , 8", s.s.


TOWEL CLAMPS AND SPONGE FORCEPS




SGI-301 Backhaus Towel Clamp 9cm 3 1/2", stainless
SGI-302 Backhaus Towel Clamp 13.5 cm 5 1/4", stainless




SGI-303 Backhaus -Roeder Towel Clamp 13.5 cm 5 1/4" with ball spots, stainless .



SGI-304 Edna Non- Perforating Towel Clamp 9.5 cm , 3 3/4", stainless
SGI-304 Edna Non- Perforating Towel Clamp 13.5 cm , 5 1/4", stainless



SGI-306 Jonest Towel Clamp
9 cm, 3 1/2", stainless




- SGI-307 Foerster Sponge Forceps
straight, smooth jaw
18 cm 7", stainless
- SGI-308 Foerster Sponge Forceps
straight, smooth jaw
24 cm 9 1/2", stainless
- SGI-309 Foerster Sponge Forceps
curved, smooth jaw
18 cm 7", stainless
- SGI-310 Foerster Sponge Forceps
curved, smooth jaw
24 cm 9 1/2", stainless
- SGI-311 Foerster Sponge Forceps
straight, serrated, 18 cm 7", s.s.
- SGI-312 Foerster Sponge Forceps
straight, serrated, 24 cm 9 1/2", s.s.
- SGI-313 Foerster Sponge Forceps
curved, serrated, 18 cm 7", s.s.
- SGI-314 Foerster Sponge Forceps
curved, serrated, 24 cm, 9 1/2", s.s.

GYNECOLOGICAL INSTRUMENT




SGI-315 Sims Uterine Sound
33 cm, 13", stainless




SGI-316 Sims Uterine Sound, malleable
33 cm, 13", silver plated




SGI-317 Hegar Uterine Dilators
Set of 13, 1/2- 25/26mm




SGI-318
Bozeman uterine Dressing forceps
straight, 25.5 cm, 10", stainless




SGI-319
Bozeman uterine Dressing forceps
curved 25.5 cm, 10", stainless




SGI-320 Schroeder Tenaculum Forceps
straight, 24 cm, 9 1/2", stainless



SGI-321 Schroeder tenaculum Forceps
curved, 25 cm, 9 1/2", stainless




SGI-322 Lahey Vulsellum forceps
16 cm 6 1/2", stainless




SGI-323
set of 8


SGI-324
set of 8



EAR AND INSTRUMENTS




- SGI-326 Ear Syringe with shield and
two tips, 2 oz (57 g) C. P.
- SGI-327 Ear Syringe with shield and
two tips, 3 oz (85 g) C. P.
- SGI-328 Ear Syringe with shield and
two tips, 4 oz (113 g) C. P.
- SGI-329 Ear Syringe with shield and
two tips, 6 oz.
- SGI-330 Ear Syringe with shield and
two tips, 8 oz.



SGI-331 Wilde Ear Forceps
12.5 cm, 5", stainless




SGI-332 Lucae Ear Forceps
12.5 cm, 5", stainless




SGI-333
Hartam Alligator Forceps
shaft 14 cm, 5 1/2"/16.5cm, 6 1/2"
serr. jaw, s.s.



SGI-334 Nasal Speculum
adult size, stainless




SGI-335 Nasal Speculum
child size, stainless




- SGI-336 Killian Nasal Speculum
18 cm, 1 3/8".
- SGI-337 Killian Nasal Speculum
5 cm, 2".
- SGI-338 Killian Nasal Speculum
7.5 cm, 3".
- SGI-339 Killian Nasal Speculum
9 cm, 3 1/2".

RETRACTORS



SGI-340 Alms Self-Retaining Retractor
7 cm, 2 3/4", stainless

SGI-341 Alms Self-Retaining Retractor
10 cm, 4", stainless




SGI-342 Weitlaner Self-Retaining Retractor
3x4 prongs, 13 cm, 5"
sharp teeth, stainless

SGI-343 Weitlaner Self-Retaining Retractor
3x4 prongs, 13 cm, 5"
blunt teeth, stainless


SGI-344 Weitlaner Self-Retaining Retractor
3x4 prongs, 16.5 cm, 6 1/2"
sharp teeth, stainless

SGI-345 Weitlaner Self-Retaining Retractor
3x4 prongs, 16.5 cm, 6 1/2"
blunt teeth, stainless




SGI-346 Senn-Muller Retractor
3 sharp prongs, 16 cm, 6 1/2", s.s.

SGI-347 Senn-Muller Retractor
3 blunt teeth prongs, 16 cm, 6 1/2", s.s.




SGI-348 Desmarres Eye Lid Retractor
14 cm, 5 1/2", size 1-4



SGI-349 Doyen-jansen Mouth Gag
10 cm, 4", stainless

SGI-350 Doyen-jansen Mouth Gag
12.5 cm, 5", stainless



SGI-351
40x40 mm

SGI-352
40x50 mm


SGI-353
50x60 mm

SGI-354
70x70 mm


- Sims
SGI-355
65x25mm
- SGI-356
70x30mm
- SGI-357
70x30mm
- 70x35m
- 85x35mm
- 90x40mm

Israel
25cm /10"


BONE CUTTERS




- SGI-358 Liston Bone Cutting Forceps
Box lock, straight 14cm, 5 1/2"s.s.
SGI-359 Liston Bone Cutting Forceps
Box lock, straight 17cm, 6 3/4"s.s.
SGI-360 Liston Bone Cutting Forceps
Box lock, straight 19cm, 7 1/2"s.s.



- SGI-361 Ruskin-Liston Bone Forceps
18.5 cm ,7 1/4",s.s.



- SGI-362 Ruskin-Liston Bone Forceps,angled,
18.5 cm ,7 1/4",s.s.



- SGI-363 Luer Bone Rongeur
straight, 18cm, 7", s.s.


PROBES AND DIRECTORS




- SGI-364 Probes D/D
12.5cm 5", stainless
SGI-365 Probes D/E
14cm 5 1/2", stainless
SGI-366 Probes D/E
15cm 6", stainless
SGI-367 Probes D/E
18cm 7", stainless
SGI-368 Probes D/E
20 cm 8", stainless
SGI-369 Probes D/E
25cm10", stainless
SGI-369 Probes D/E
25cm10", stainless



- SGI-370 Probes with Eye
14cm 5 1/2", stainless
SGI-371 Probes with Eye
15cm 6", stainless
SGI-372 Probes with Eye
20cm 8", stainless
SGI-373 Probes with Eye
25cm 10", stainless




- SGI-374 Williams Lacrimal Probes, size 0/00.
SGI-375 Williams Lacrimal Probes, size 1/2.
SGI-376 Williams Lacrimal Probes, size 3/4.
SGI-377 Williams Lacrimal Probes, size5/6.
SGI-378 Williams Lacrimal Probes, size6/7.



- SGI-379 Grooved Director
14 cm, 5 1/2,stainless



- SGI-380 Grooved Director probe point
14 cm, 5 1/2,stainless




- SGI-381 Mall Probe and Seeker
15 cm, 6,stainless


SCALPEL HANDLES



- SGI-382 Scalpel Handles
#3, 12.5 cm,5", stainless




- SGI-383 Scalpel Handles
#4, 14 cm, 5 1/2", stainless



- SGI-384 Scalpel Handles
#7k, 16 cm, 6 1/4", stainless



- SGI-385 Scalpel Handles
#9, 12.5 cm,5", stainless



- SGI-386 Scalpel Handles
#3L, 20 cm,8", stainless



- SGI-387 Scalpel Handles
#4L, 22 cm,8 1/2", stainless

DENTAL INSTRUMENTS




- SGI-388 Dental Tweezer
14cm, 5 1/2, stainless
SGI-389 Dental Tweezer
15cm, 6", stainless



- SGI-390 Dental Tweezer
14cm, 5 1/2", stainless
SGI-391 Dental Tweezer
15cm, 6", stainless




- SGI-392 Tooth Extracotor, stainless



- SGI-393 Tooth Extracotor, stainless



- SGI-394 Tooth Extracotor, stainless




- SGI-395 Tooth Extracotor, stainless




- SGI-396 Root Extracotor, stainless


CHIROPODY INSTRUMENT




- SGI-397 Nail Nipper with double pring
14 cm , 5 1/2", stainless




- SGI-398 Nail Nipper with barrel pring
12 cm , 4 1/2", chrome
SGI-399 Nail Nipper with barrel pring
12.5 cm ,5", chrome




- SGI-400 Ingrown Toenail splitter
with spring, 12 cm , 4 1/2",chrome




SGI-401 Nail clipper straight edge, 14cm , 5 1/2, s.s.




SGI-402 Nail clipper Scissors Type 10.5 cm, 4 1/2".



SGI-403 Nail Nipper Angled Blade 11.5 cm , 4 1/2".



SGI-404 Nail Nipper Curved Blade 14 cm , 5 1/2".




SGI-405 Toenail Clipper
SGI-406 Finger Nail Clipper.




SGI-407 Nail file, 14cm, 5 1/2", s.s.




SGI-408 Corn Cutter With Blade




SGI-409 Nail Scissors straight 9 cm, 3 1/2, s.s.




SGI-410 Nail Scissors curved 9 cm, 3 1/2, s.s.




SGI-411 Nail Scissors straight, 9 cm, 3 1/2",stainless




SGI-412 Nail Scissors fine, 9 cm, 3 1/2",stainle



SGI-413 Nail Scissors 9 cm, 3 1/2",stainle



SGI-414 Cuticle Scissors 9 cm, 3 1/2",chrome
SGI-415 Cuticle Scissors 9 cm, 3 1/2",stainless



SGI-416 Nail Scissors 7.5 cm, 3",stainless




SGI-417 Littauer Cilla Forceps 9cm 3 1/2", 4mm wide jaws, fine serrations.




SGI-418 Bergh Cilla Forceps 9cm 3 1/2", 5mm wide angular jaws, fine serrations.

TUBE OCCLUDING FORCEPS AND SUCTION TUBES




SGI-419 Tube Occluding Forceps angled smooth jaw, 11.5 cm 4 1/2",s.s.



SGI-420 Tube Occluding Forceps straight smooth jaw, 11.5 cm 4 1/2",s.s.


SGI-421 Tube Occluding Forceps straight smooth jaw, 18 cm 7",s.s.




SGI-422 Tube Occluding Forceps 7 1/2 narrow jaw, straight ,serrated, stainless. lock guard prevents damage to tubing by lock joint.




SGI-423 Poole Suction Tube , straight 30 fr, 22cm 8 1/2",s.s.



SGI-424 Poole Suction Tube , curved 23 fr, 22cm 8 1/2",s.s.




SGI-425 Yankauer Suction Tube with threaded end



SGI-426 Yankauer Suction Tube with slip-on end

BONE SAWS




SGI-427 Satterlee Bone Saw,31cm




SGI-428 Langenback Bone Saw,31cm


MISCELLANEOUS INSTRUMENTS




SGI-429 Ward Scissors, 18 cm , 7", c.p.
SGI-430 Ward Scissors, 20 cm , 8", c.p.




SGI-439 Northbent Stitch Scissors, curved
12cm , 4 3/4",s/s




SGI-431 smith Scissors, 18 cm , 7", c.p.
SGI-432 smith Scissors, 20 cm , 8", c.p.




SGI-440 O,Brien Stitch Scissors,
delicate, angled
9cm, 3 3/4",s/s




SGI-433 Straight Trimmers,
18 cm , 7", c.p.




Cusco
SGI -441(75x32mm)
SGI-442(85x35mm)
SGI-443(100x37mm)




SGI-434 Bent Trimmers,
20 cm , 8", s.s.




SGI-444 Graves Vaginal Speculum
Large, 11.5x3.5cm blade ,S/S
SGI-445 Graves Vaginal Speculum
Medium, 10x3cm blade ,S/S
SGI-446 Graves Vaginal Speculum
Small, 7.5x2 cm blade ,S/S




SGI-435 Housedhold Scissors,
15 cm , 6", s.s.



SGI-437 Probe Scissors,13 cm , 5", s.s.




SGI-447 Pederson Vaginal Speculum
Large, 12x2.5cm blade ,S/S
SGI-448 Pederson Vaginal Speculum
Medium, 10x2.2cm blade ,S/S
SGI-449 Pederson Vaginal Speculum
Small, 7.5x1.3cm blade ,S/S



SGI-438 Shortbent Stitch Scissors, curved
8.5cm , 3 1/2",s/s

SURGICAL BLADES




10 11 12 15




20 21 22 23

- SGI-450 No. 10 Surgical Blades c.s.
- SGI-451 No. 11 Surgical Blades c.s.
- SGI-452 No. 12 Surgical Blades c.s.
- SGI-453 No. 15 Surgical Blades c.s.
- SGI-454 No. 20 Surgical Blades c.s.
- SGI-455 No. 21 Surgical Blades c.s.
- SGI-456 No. 22 Surgical Blades c.s.
- SGI-457 No. 23 Surgical Blades c.s.

TUBE GAUZE APPLICATORS



SGI-458 Gauze Applicator size C(01)



SGI-459 Gauze Applicator size D(02)

Salim



Surgical Instruments



SALIM

Group Of Industries

P.O.Box: 2561, Sialkot -Pakistan.

Tel : 0092-52-4597649

Fax : 0092-52-4587056

E-Mails: salimgroup@brain.net.pk

: malikshahidsaleem@hotmail.com

Web : www.salim-group.com

malikshahidsaleem@hotmail.com